

Competency Checklist



This course will use the *epoc Blood Analysis System with NXS Host Resource Guide* and *PEP Connect* as training resources.

The student places a checkmark beside the competency when it is completed. When all competencies are checked, the trainer and the student sign and date below as a record of completion.

Topic/Tasks	Completed
System Overview	
Identify major hardware components of the epoc System	
Identify major software functions of the epoc System	
Sample Processing	
Use the NXS Host to login to the epoc Application	
Use the Reader to insert a Test Card	
Use the NXS Host to enter test information	
Use a syringe to inject a sample	
Use a Care-Fill™ capillary tube to inject a sample	
Review test results	
Quality Assurance	
Review Electronic Quality Control (EQC)	
Troubleshooting	
Discuss troubleshooting common result errors	
Discuss troubleshooting common hardware errors	

Competency Checklist



Trainer:

Student:

Date:

Trainer comments