

# Dimension Vista® Systems Virtual Training

## Preparation Checklist



### General information about the training:

- Participants should receive a Registration Confirmation email for the upcoming virtual training class that includes the following:
  - Weblink to join the virtual training session
  - Link to the Course Pre-Training (*if applicable*)
  - Link to download training materials
- Registered participants must be present for the entire duration of the course to receive a certificate of completion.
- Participants will require the use of a computer with speaker and microphone. In the event that the computer does not have an audio option, we will provide a dial in option at the beginning of the session.
- An Installed and Validated Instrument is required for the duration of the virtual training event.

Training Materials Received from Siemens Healthineers	Completed
Dimension Vista Systems Quick Reference Guide	
Dimension Vista® Systems Virtual Training Workbook	
Orange folder containing various documents	
Notepad	
Pen	
<b>Reagent/QC Needs</b>	
GLUC, BUN, URCA flexes	
Daily Quality Control	
CHEM I Calibrator with lot specific package insert (in box)	
<b>Consumables Needs</b>	
Aliquot plates, Cuvettes Cartridge, LOCI Reaction Vessels	
Bulk Fluids as needed	
<b>Samples</b>	
Dimension Sample Cups, Sample tubes, SSC's	
Vista Sample Racks	
5-10 patient samples, barcodes if available	
<b>Maintenance/Troubleshooting Needs</b>	
Cleaning Supplies: Gauze and Distilled Water, alcohol prep pad	
Flex inserts- R1 & R3	
Air Filters, Drain brush	
WPM: Bleach Tablets	
Wrenches: 7 & 8 mm	