

# ADVIA® 2120i Hematology System Virtual Training

## Preparation Checklist



### General information about the training:

- Participants should receive a Registration Confirmation email for the upcoming virtual training class that includes the following:
  - Weblink to join the virtual training session
  - Link to the Course Pre-Training (*if applicable*)
  - Link to download training materials
- Registered participants must be present for the entire duration of the course to receive a certificate of completion.
- Participants will require the use of a computer with speaker and microphone. In the event that the computer does not have an audio option, we will provide a dial in option at the beginning of the session.
- An Installed and Validated Instrument is required for the duration of the virtual training event.

Training Materials Received from Siemens Healthineers	Completed
ADVIA 2120i Hematology System Quick Reference Guide	
ADVIA 2120i Hematology System Virtual Training Workbook	
Orange folder containing various documents	
Notepad	
Pen	
Consumables/Supplies Received from Siemens Healthineers	
ADVIA 2120/2120i RBC Flow Cell Wash (SMN: 10734561)	
ADVIA 2120/2120i Perox Flow Cell Wash (SMN: 11306489)	
ADVIA 2120/2120i Aspiration Pathway Wash (SMN: 11306490)	
ADVIA 2120/2120i Vent Line Wash (SMN: 11306492)	
ADVIA OPTIpoint™ (SMN:10312283)	
Tubing Standard 0.81-vent line (SMN: 10310236)	
Quality Control Needs	
ADVIA TESTpoint™ 3-in1, ADVIA TESTpoint™, or some other hematology QC product	
Sample Processing Needs	
6 barcoded purple tops with whole blood	
2 non-barcoded purple tops with whole blood *if no blood samples can be acquired, may need to use QC blood to make samples	
How the System Works Needs	
1-clot filter (SMN: 10310090)	

<b>Scheduled Maintenance Needs</b>	
DI water squirt bottle	
Lint free tissues	
Beaker of 6% bleach (sodium hypochlorite)	
Beaker of 25% solution of 6% bleach and DI water	
Beaker of EZ Wash	
Red needle cover	
Cotton swabs	
Stylet	
Flowcell cleaning syringe (SMN: 10316451)	
TUBING TYGON .030 ID RED (SMN: 10310229)	
Parker Super O-Lube	
RBC 0.22mm Sheath Filter (SMN: 10311257)	
Perox 0.45mm Sheath Filter (SMN: 10488816)	
1 box of ADVIA 2120/2120i RBC Flow Cell Wash (SMN: 10734561)	
1 box of ADVIA 2120/2120i Perox Flow Cell Wash (SMN: 11306489)	
1 box of ADVIA 2120/2120i Aspiration Pathway Wash (SMN: 11306490)	
1 box of ADVIA 2120/2120i Vent Line Wash (SMN: 11306492)	
0.081 Vent line tubing- (SMN:10310236)- need two separate 12 inch pieces	
<b>Gains Needs</b>	
ADVIA Optipoint	