

Viva-ProE® System Virtual Training

Preparation Checklist



General information about the training:

- Participants should receive a Registration Confirmation email for the upcoming virtual training class that includes the following:
 - Weblink to join the virtual training session
 - Link to the Course Pre-Training (*if applicable*)
 - Link to download training materials
- Registered participants must be present for the entire duration of the course to receive a certificate of completion.
- Participants will require the use of a computer with speaker and microphone. In the event that the computer does not have an audio option, we will provide a dial in option at the beginning of the session.
- An Installed and Validated Instrument is required for the duration of the virtual training event.

Training Materials Received from Siemens Healthineers	Completed
Viva-ProE System Quick Reference Guide	
Viva-ProE System Virtual Training Workbook	
EMIT Drugs of Abuse Cross-Reactivity List	
Orange folder containing various documents	
Notepad	
Pen	
Reagents/System Liquids	
Tests bottles on reagent rotor	
Calibrators and Controls for tests on reagent rotor	
Sodium Hypochlorite Solution (cleaning solution)	
Hydrochloric Acid Solution (0.1N HCl)	
System Solution + graduated cylinder to measure	
Coolant fluid + graduated cylinder to measure	
Consumables Needs	
2 mL cups	
Sample tubes	
Transfer pipettes	
Metal adaptors	
Sample tube rack	
Samples	
5 urine samples	

Maintenance/Troubleshooting Needs	
USB for Backups	
Maintenance Kit	
Small Phillips head screwdriver + graduated cylinder to measure	